

HITACHI 911

KREATININ ENZIMATIKUS stabil folyékony
(46161-46162)

ST	[KREE]					
ASSAY CODE	[2poin End] [10] []		WAVELENGTH (2nd/Primary):	[700]/[546]		
ASSAY POINT	[16] - [34] - [0] - [0]		Diluent/Rgt Stability	[*] [99]		
S. VOL (Regular)	[8] [0]	[0]				* Diluent number
S. VOL (Decrease)	[4] [0]	[0]				
S. VOL (Increase)	[16] [0]	[0]				
ABSORB LIM	[-32000]			[32000] [Increase]		
PROZONE LIM	[-32000]			[32000] [Upper]		
REAGENT	R1 [180] [0]	[**]	[0]			
	R2 [] [0]	[]	[0]			
	R3 [60] [0]	[**]	[0]			
	R4 [] [0]	[]	[0]			
CALIB TYPE	[Linear]		[2]-[2] [0] []			
AUTO TIME OUT BLANK	[0]		SD LIM	[999]		
	SPAN [0]		DUPLICATE LIM	[32000]		
	2POINT [0]		SENSITIVITY LIM	[0]		
	FULL [0]		S1ABS LIM	[-32000] [32000]		
AUTO CHANGE OF LOT	[cancel]		COMPENSATED LIM	[]		
CHANGE OF BOTTLE	[cancel]					
TEST	[KREE] [**]		TEST NAME	[KREAE]		UNIT [umol/l]
DATA MODE	[On Board]		REPORT NAME	[KREAE]		
QC INTERVAL	[0]		INSTRUMENT FACTOR (y=ax +b)	a[1.0] b[0.0]		
EXPECTED VALUES	[40] - [100] umol/l					
TECHNICAL LIMIT	[5] - [1770]					
STD	CONC.	POS	SAMPLE	PRE	DIL	CALIB
(1)	[0.0]	***	[10]	[0]	[0]	[**]
(2)	[***]	***	[10]	[0]	[0]	[**]
(3)	[***]	***	[10]	[0]	[0]	[**]
(4)	[***]	***	[10]	[0]	[0]	[**]
(5)	[***]	***	[10]	[0]	[0]	[**]