



**Adenovirus Detection Kit v1**  
**USER MANUAL**

For *in vitro* Diagnostic Use



Document Code: MB83v1f  
Approval Date: September 2013

IVD



## Contents

	<u>Page</u>
1. Product Description	1
2. Content	1
3. Storage	1
4. Required Materials and Devices	1
5. Important Notes and Safety Instructions	2
6. Product Use Limitations	2
7. Pathogen	2
8. Method	3
9. Procedure	4
9.1. DNA Isolation	4
9.2. Kit Components	4
9.2.1. PCR Master Mix	4
9.2.2. Internal Control	4
9.2.3. Positive Control	5
9.3. Preparing the PCR	5
9.4. Programming the Real-Time PCR Instrument	5
10. Analysis	6
11. Specifications	7
11.1. Sensitivity	7
11.2. Cross-Reactivity	7
11.3. Reproducibility	7
12. References	8
13. Symbols	8
14. Contact Information	9



## 1. PRODUCT DESCRIPTION

Bosphore® Adenovirus Detection Kit v1 detects Adenovirus DNA in human serum or plasma, encompassing all the major Adenovirus subtypes. The analytic sensitivity is 88 copies/rxn. A region within the hexon gene is amplified and fluorescence detection is accomplished using the FAM filter.

An internal control has been integrated into the kit in order to check PCR inhibition. The amplification data of the internal control is detected with the Cy5/HEX/JOE filter. The internal control can be added either during DNA extraction or PCR step.

## 2. CONTENT

Bosphore® Adenovirus Detection Kit v1 is composed of Real-Time PCR reagents:

Component	REAGENT	100 Reactions	50 Reactions	25 Reactions
1	dH <sub>2</sub> O	(1000 µl)	(500 µl)	(500 µl)
2	PCR Master Mix	(1628 µl)	(814 µl)	(407 µl)
3	Internal Control	(22 µl)	(15 µl)	(15 µl)
4	Positive Control 1	(88µl)	(44 µl)	(22 µl)

## 3. STORAGE

Bosphore® Adenovirus Detection Kit v1 PCR reagents should be stored at -20°C. Repeated thawing and freezing (>3x) should be avoided since it may reduce sensitivity. If the components are to be used in small amounts, they should be frozen in aliquots.

While preparing the PCR; the components should not be exposed to room temperature for more than 10 min. and the detection mix components should not be exposed to light or air more than necessary, vials must be kept closed except during pipetting. We recommend preparing the PCR on a cooling block and keeping the detection mixes within a closed container.

The components maintain their stability until the expiry dates on the labels, if they are stored at advised conditions.

## 4. REQUIRED MATERIALS AND DEVICES

- Montania® 483 or Montania® 4896 Real-Time PCR Instrument (Anatolia Geneworks), or another Real-Time PCR system with FAM and Cy5/HEX/JOE filters (such as iCycler, iQ5, CFX96–BioRad, LightCycler 2.0, 480-Roche, 7300, 7500 Real-Time PCR System-ABI, Stratagene Mx3005P, Mx3000P-Agilent, LineGeneK, LineGene 9600-Bioer, Rotorgene 6000, Q-Qiagen)
- 0.2 ml Thin-Wall PCR tubes, PCR plates or strips
- Magnesia® 16 Nucleic Acid Extraction System / Magnesia® Viral Nucleic Acid Extraction Kit / Bosphore® Viral DNA Extraction Spin Kit (Anatolia Geneworks) or other high quality viral DNA extraction kits and systems
- Deep freezer (-20°C)

- Desktop centrifuge with rotor for 2 ml. microcentrifuge tubes
- Calibrated adjustable micropipettes
- DNase, RNase, pyrogen free micropipette tips with filters
- DNase, RNase, pyrogen free 1.5 or 2 ml. microcentrifuge tubes
- Disposable laboratory gloves

## 5. IMPORTANT NOTES AND SAFETY INSTRUCTIONS

### Important!:

- The product should be delivered on dry ice. Check for presence of dry ice upon arrival.
- Check for the expiry dates on the box and tube labels, upon arrival. Do not use expired products or components.
- Calibrated or verified micropipettes, DNase, RNase, pyrogen free micropipette tips with filters, and DNase, RNase, pyrogen free microcentrifuge tubes should be used.
- Before starting a test procedure, all components should be thoroughly thawed. After thawing, all components should be centrifuged briefly (spin-down for 3-5 seconds), and mixed well to ensure homogeneity prior to use.
- The kit components should be kept on ice or a cooling block until the reaction is prepared, and they should be quickly returned to -20°C.
- PCR and nucleic acid isolation must be performed in different compartments. Samples should be stored separately to avoid contact with the kit components.
- Pathogen information should be reviewed to be aware of the health related risks.
- Serum samples should be handled with extreme caution, suitable class microbiological safety cabinet should be used: Physical contact with pathogens should be avoided by; wearing lab coats and gloves, no allowance for eating or drinking within the workspace, prevention of unauthorized individuals' access to the working area.
- All the pathogenic wastes produced during the nucleic acid isolation step; including the serum samples and material contacted with them, should be discarded into medical waste and disposed safely.

## 6. PRODUCT USE LIMITATIONS

- All the components may exclusively be used for in vitro diagnostics.
- This product should be used in accordance with this user manual, by personnel specially trained to perform in vitro diagnostic procedures.

## 7. PATHOGEN

### Causative Agents

Adenovirus is a relatively large non-enveloped dsDNA virus which is a member of the Adenoviridae family. They are grouped into seven species (A to G) and 53 serotypes based on immunochemical responses,

nucleic acid characteristics, hexon and fiber protein characteristics, biological properties and phylogenetic relationships. The virus is able to replicate in the nucleus of mammalian cells using the host's replication machinery. It mostly results in infections of the upper respiratory tract but AdV types 40 and 41 can also cause gastroenteritis. (1,2)

### **Epidemiology**

Human adenoviruses (HAdVs) were the first respiratory viruses to be isolated and characterized. The viruses are associated with acute respiratory disease (ARD), conjunctivitis, genitourinary infections, and gastroenteritis. Adenovirus is a very common infection and responsible for between 2% and 5% of all respiratory infections. It usually occurs during childhood but all ages are susceptible. In winter, infection with type 4 or 7 causes recognizable illness in military recruits, with about 25% requiring hospitalization for fever and lower respiratory tract disease (3, 4)

### **Modes of Transmission:**

AdV are usually transmitted via the fecal-oral route and the respiratory route. Close contact with infected people such as touching or shaking hands, inhaling the air by coughing and sneezing, touching an object or surface with adenoviruses on it then touching your mouth, nose, or eyes before washing your hands can cause the viral transmission. (4)

## **8. METHOD**

Bosphore® Adenovirus Detection Kit v1 is based on the Real-Time PCR method. Polymerase chain reaction is a technique that is used for amplification of a DNA region. The reaction occurs by the repeating cycles of heating and cooling. The main components of PCR are primers, dNTPs, Taq polymerase enzyme, buffer solution and template. As a brief explanation, primers are small synthetic DNA those anneal to the specific regions of the template in order to start the synthesis. dNTPs are the building blocks of the amplified products. Taq polymerase amplifies the DNA template. Buffer solution provides the pH adjustment required for the reaction and template, as referred, is the target region for synthesis.

In Real Time PCR technique, in contrast to conventional PCR, PCR product can be monitored during the reaction. Therefore Real-Time PCR obviates the need for further analysis methods like gel electrophoresis, whereby minimizing the risk of contamination. Dual labeled probes employed in the reaction in addition to the conventional PCR reagents, enable detection of the amplified target with increased sensitivity.

The assay utilizes the 5' exonuclease activity of Taq Polymerase to cleave a dual-labeled fluorescent probe during the extension phase of PCR. The probe is labeled at the 5' end with a fluorescent 'reporter' molecule, and at the 3' end with another fluorescent molecule that acts as a 'quencher' for the 'reporter'. When the two fluorophores are in close proximity, and the reporter is excited by light, no reporter fluorescence can be detected. During the elongation step of PCR, Taq Polymerase encounters and cleaves the probe bound to the template. As the reporter is freed from the suppressing effect of the quencher, fluorescence signal can be detected.

The fluorescence generated by the reporter increases as the PCR product is accumulated; the point at which the signal rises above background level and becomes distinguishable, is called the threshold cycle ( $C_T$ ). There is a linear relationship between the log of the starting amount of a template and its threshold cycle, thus starting amount of unknown templates can be determined using standard curves constructed using  $C_T$  values of the known starting amounts of target templates.

Bosphore® Adenovirus Detection Kit v1 employs multiplex PCR, and an internal control is incorporated into the system in order to control the isolation procedure and to check for possible PCR inhibition. Adenovirus DNA and an internal control are co-amplified in a single reaction, using sequence-specific primers. The fluorescent signal generated by the Adenovirus amplification is detected by a probe labeled at the 3' end with FAM, through the FAM channel. The fluorescent signal generated by the internal control amplification, is detected by a second probe (labeled at the 5' end with a different reporter molecule, Cy5/HEX/JOE) through the Cy5/HEX/JOE channel.

## **9. PROCEDURE**

### **9.1. DNA Isolation**

We recommend that the Magnesia® 16 Nucleic Acid Extraction System / Magnesia® Viral Nucleic Acid Extraction Kit / Bosphore® Viral DNA Extraction Spin Kit (Anatolia Geneworks) or another high quality isolation system is used with Bosphore® Adenovirus Detection Kit v1. The DNA isolation should be performed according to the manufacturers' instructions. The starting volume for Anatolia Geneworks systems is 400 µl, the elution volume is 60 µl and the amount of internal control that should be used during isolation for each system is 5 µl.

### **9.2. Kit Components**

#### **9.2.1. PCR Master Mix**

HotStarTaq DNA Polymerase: HotStarTaq DNA Polymerase is a modified form of a recombinant 94 kDa DNA polymerase, originally isolated from *Thermus aquaticus*, cloned into E.Coli. The enzyme is provided in an inactive form. It is activated by a 15-minute 95 °C incubation step. This prevents the formation of misprimed products and primer-dimers during reaction setup and the first denaturation step, leading to high PCR specificity and accurate Detection. PCR Buffer: contains Tris-Cl, KCl,  $(NH_4)_2SO_4$ , 8 mM  $MgCl_2$ , pH 8.7 (20°C). dNTP Mix: Contains ultrapure quality dATP, dGTP, dCTP ve dTTP/dUTP.

PCR Master Mix also contains Adenovirus-specific forward and reverse primers and a dual-labeled probe, and the internal control-specific forward and reverse primers and a dual-labeled probe.

#### **9.2.2. Internal Control**

An internal control is included in the kit to control PCR inhibition. The internal control is a synthetic DNA molecule. It is added directly into the PCR master mix to control the PCR inhibition. For this purpose, 0.2 µl of internal control should be added for each reaction into the master mix. Lack of internal control amplification in the FAM negative samples, may indicate PCR inhibition. In this case, isolation and PCR

should be repeated. In samples that contain a high viral load, internal control can be suppressed and no increase of the signal is detected. Please use the table below for the interpretation of internal control data:

Adenovirus (FAM)	Internal Control (Cy5/HEX/JOE)	Interpretation
+	+	Sample positive
-	+	Sample negative
+	-	Sample positive
-	-	Repeat the test!

### 9.2.3. Positive Control

The positive control contains Adenovirus DNA. It should be included in the PCR to test the efficiency of the PCR exclusively. The threshold cycle for the positive control is given in the acceptance criteria table (Section 10. Analysis). Threshold cycles higher than the acceptance criteria may indicate an efficiency loss in the reaction.

### 9.3. Preparing the PCR

The positive and the negative control (Extracted Adenovirus negative serum/plasma) should be added into the PCR reaction together with the samples. Make sure that all the kit components are thawed before use. Refer to the table below for preparing the PCR. It is for only one reaction, multiply these values with the sample number to find the values required for the master mix. While preparing master mixes for more than 5 samples, an extra 10% should be added to the total sample number.

<b>PCR Master Mix</b>	14.8 µl
<b>Internal Control</b>	0.2 µl
<b>Sample DNA (Negative/Positive Control)</b>	10 µl
<b>Total Volume</b>	25 µl

Pipette 15 µl of the master mix into the PCR tubes or strips, and add 10 µl of DNA (sample/positive or negative control). Close the tube cap. Make sure that the solution in each tube is at the bottom of the tube. Centrifuge if necessary.

### 9.4. Programming the Real-Time PCR Instrument

The thermal protocol for Bosphore® Adenovirus Detection Kit v1 is composed of an initial denaturation for activation the HotStarTaq DNA Polymerase, a two-step amplification cycle and a terminal hold. The real-time data is collected at the second step of the amplification cycle.

Initial denaturation	95°C	14:30 min.	} 50 cycles
Denaturation	97°C	00:30 min.	
Annealing (Data Collection)	53°C	01:30 min.	
Synthesis	72°C	0:15 min.	
Hold	22°C	05:00 min.	



To start a Real-Time PCR reaction using the Bosphore® Kits, the following steps should be completed:

- Choose the filter pairs to be used (FAM and Cy5/HEX/JOE),
- Identify unknown samples, positive and negative controls, Select the correct thermal protocol.
- Start the protocol

## 10. ANALYSIS

By the end of the thermal protocol, the Real-Time PCR Instrument software automatically calculates the baseline cycles and the threshold.

Example of an amplification curve is given in Fig. 1.




Fig. 1: Amplification Curve of a Bosphore® Adenovirus test

Analysis of the results should be performed by trained personnel who have received the required training for analyzing Real-Time PCR data. We recommend that the test results must be evaluated by an expert clinician, taking the patient's clinical findings and the results of other tests into consideration.

All analysis is done automatically in routine use. However, when the trained personnel, who have received the required training from manufacturer, consider it as necessary, and if the system's software allows pulling down the threshold as much as possible in order to detect low positive samples, attention should be paid to keep the threshold line above the background and to keep the correlation coefficient at the maximum possible value (and within its acceptance criteria).

Positive control of Bosphore® Adenovirus v1 is essential for accurate result analysis. The cycle threshold acceptance criterion for the positive control is  $32 \pm 4$ . Test results should not be reported unless the assay results meet the criteria stated above. Please contact the manufacturer if impairment in the product's performance is observed (See the last page for contact information).

The qualitative results of the test are displayed on the "Report Mode" screen. The samples that cross the threshold in FAM channel are displayed as positive whereas samples that do not cut the threshold are

displayed as "Negative" or "No Ct". These samples are regarded as negative or having a bacterial load below the detection limit of the assay. For these undetectable samples, the CY5/HEX/JOE data of the internal control should also be checked to avoid false negative results.

The following table shows the possible results and their interpretation:

Signal detected in FAM filter pair	The sample contains Adenovirus DNA, the result is positive	No need to check the internal control since the sample is positive (high positive samples may suppress the signal from the internal control)
No signal in FAM, signal in Cy5/HEX/JOE	The Adenovirus DNA in the sample is not detectable	Signal from Cy5/HEX/JOE filter pair rules out the possibility of PCR inhibition
No signal in FAM and Cy5/HEX/JOE	The diagnosis is inconclusive	No signal in Cy5/HEX/JOE points out to PCR inhibition or to a problem in DNA isolation

## 11. SPECIFICATIONS

### 11.1. Sensitivity

Analytical sensitivity may be expressed as the limit of detection: i.e. the smallest amount of the target marker that can be precisely detected. The detection limit of an individual analytical procedure is the lowest amount of nucleic acid in a sample which can be detected but not necessarily quantitated as an exact value. The analytical sensitivity or detection limit for NAT assays is expressed by the 95% positive cut-off value.

The analytical detection limit for Bosphore® Adenovirus v1 was found to be 88 copies/rxn (p=0.05). The sensitivity was determined using serial dilutions of DNA calibrated with the previously quantitated Adenovirus DNA Control. The dilutions were tested in different runs in replicates. The results were analyzed by probit method.

### 11.2. Cross-Reactivity

To eliminate potential cross-reactivity, both assay design evidence and experimental studies were employed. Primer and probe sequences were checked for possible homology to other known pathogen sequences by sequence comparison analysis using database alignment. Samples of CMV, EBV, MTBC, Parvovirus B19, BKV with known high positivity were tested, and found negative.

### 11.3 Reproducibility

Reproducibility data (on  $C_T$  value basis) were obtained by the analysis of positive control of the Bosphore® Adenovirus Detection Kit v1. Test was performed in at least 4 replicates by 3 different operators, on multiple days, using 3 different lots. The resulting data is given in Table 1:

Table 1: Reproducibility Data.

<b>ADV (7000 cop/rxn)</b>	<b>Standard deviation</b>	<b>Variance</b>	<b>Coefficient of variation [%]</b>
<b>Intra-assay Variability N=4</b>	0,051	0,003	0,160
<b>Inter-lot Variability N=3</b>	0,071	0,005	0,222
<b>Inter-operator Variability N=3</b>	0,056	0,003	0,175
<b>Total Inter-assay Variability N=5</b>	0,053	0,003	0,166

## 12. REFERENCES

1. Smith J.G., Wiethoff C.M., Stewart P.L., and Nemerow G.R., Adenovirus, Curr Top Microbiol Immunol. 2010; 343: 195–224.
2. Metzgar D., Skochko G., Gibbins C., Hudson N., Lott L., Jones M.S., Evaluation and Validation of a Real-Time PCR Assay for Detection and Quantitation of Human Adenovirus 14 from Clinical Samples, 2009, PLoS ONE 4(9): e7081. doi:10.1371/journal.pone.0007081
3. Jothikumar N., Cromeans T.L., Hill V.R., Lu X., Sobsey M.D., and Erdman D.D., Quantitative Real-Time PCR Assays for Detection of Human Adenoviruses and Identification of Serotypes 40 and 41, Applied and Environmental Microbiology, June 2005, p. 3131–3136
4. Berkow, R., Beers, M.H. and Burs, M.H. (2004) Infectious Diseases. In: Merck Manual Diagnosis and Therapy. Seventeenth Edition. Merck & Co. Inc. New Jersey, USA

## 13. SYMBOLS



Use by



Lot/Batch



Catalog number



Temperature limitation



Caution, consult accompanying documents



Manufacturer



In Vitro Diagnostic Medical Device

#### 14. CONTACT INFORMATION



**Anatolia Tani ve Biyoteknoloji A.Ş.**

Egitim Mh. Kasap Ismail Sk.

No:10/23 Kadikoy 34722

ISTANBUL-TURKEY

Phone: +90 216 330 04 55

Fax: +90 216 330 00 42

E-mail: [info@anatoliagenetworks.com](mailto:info@anatoliagenetworks.com)

[www.anatoliagenetworks.com](http://www.anatoliagenetworks.com)

RegisteredTrademarks: Anatolia Geneworks® Montania®, Magnesia® and Bosphore® are registered trademarks of Anatolia Tani ve Biyoteknoloji A.S.