# MASTAZYME HSV-1/2

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	<b>Contents</b>	<b>Page</b>
1.	Intended Use	10
2.	Introduction	10
3.	Principle of the Test	11
4.	Kit Contents	11
5.	Materials Required but not Provided	12
6.	Warnings and Precautions	12
7.	Storage and Stability	12
8.	Specimen Collection and Handling	13
9.	Assay Procedure	13
10.	Results and Interpretation	14
11.	Assay Performance	15
12.	References	15

## 1. Intended Use

MASTAZYME HSV-1/2 (Herpes simplex 1/2) Antibody ELISA has been designed for the detection and the quantification of specific IgG, IgM or IgA antibodies respectively against Herpes simplex 1/2 in serum and plasma. Further applications in other body fluids are possible and can be provided on request.

This assay is intended for *in-vitro* diagnostic use only.

All laboratory test results should be interpreted in conjunction with other clinical data. The clinical judgement and further tests have to be taken into account additionally.

## 2. Introduction

Herpes simplex virus type 1 and type 2 are ubiquitous pathogens of humans that usually causes either asymptomatic infections or mild skin and mucosal diseases. Antibodies to HSV-1 occur in about 90 % of adults. Normally HSV-1 is transmitted by oral secretions or open wounds prior to the age of five. Recently in adults primary infections were observed, too.

After the primary infection some viruses establish a latent state in their host cells (mostly ganglial cells). The virus DNA is integrated into the genome of the host cell where it remains until the infected person dies. After stimulation of the host cell recurrent infection occurs which is called an exacerbation when clinical symptoms appear. The recurrence may be caused by different kinds of traumas, as fever or physiological changes and diseases. Immunosuppressed people may show a severe clinical course.

HSV-1 causes different clinical symptoms in about 10 % of the primary infections. The major clinical manifestations associated with HSV-1 infections are gingivostomatitis, keratitis, conjunctivitis, vesicular eruptions of the skin, encephalitis, eczema and some lethal infections of newborns. HSV-1 causes 85 % and HSV-2 15 % of oral primary infections. Recurrent infection occurs in form of labial fever blisters. After ulceration and scabbing of these blisters complete recovery occurs within 10 days. The central nervous system may be involved in both primary and recurrent infections. In some cases HSV-1 infection leads to a meningitis with different neurological symptoms. People at an increased risk for serious or prolonged HSV infections are those with eczema, severe burns or a defect in their cell-mediated immunity.

Antibodies to HSV-2 occur in about 20 % of all adults. In lower social classes and in sexually promiscuous people the prevalence is higher (about 60 %). HSV-2 causes different clinical symptoms. Of major importance is the Herpes genitalis syndrome which occurs principally in adults. The preceding primary infection will be transmitted via sexual contact. Typical loci of appearance are the mucosas of the human genital tract. Additionally HSV-2 is among viruses suspected of inducing cervix carcinoma in women. In some patients contiguous skin regions are involved particularly on the buttocks or perianal area. The Herpes genitalis exacerbations are normally endogenic relapses which lead to the same blistered symptoms like Herpes labialis. In some cases a HSV-2 caused meningitis occurs the course of which is much milder than a HSV-1 caused encephalitis and which symptoms are always reversible. This meningitis appears especially in connection with a HSV-2 primary infection. The most severe complication of genital HSV infection is the neonatal disease which is caused by an infection during or shortly after the delivery. The risk of neonatal infections seems to be higher during symptomatic primary infection of the mother.

The common manifestations of HSV infections are so typical that the infection may be easily diagnosed on clinical recognition alone. The "gold standard" for diagnosis of HSV infection remains isolation of the virus in tissue culture.

Diagnosis of the primary infection by HSV-1/HSV-2 can be confirmed by a significant rise of the IgG titer within 6 to 10 days. A finished infection can be monitored by the IgG ELISA. In case of a suspicion of HSV encephalopathy it is recommended to perform a parallel determination of both HSV-specific antibodies (IgG and IgM) in serum and liquor.

### 3. Principle of the Test

The principle of the test reaction can be described in four stages.

#### 3.1 Serum incubation

Specific antibodies bind to the antigens on the solid phase to form a stable immune complex. After a 60 minutes incubation at room temperature the wells are washed with prediluted wash buffer to remove all non-reactive serum components.

#### 3.2 Conjugate incubation

The anti-human-IgG /-IgM /-IgA horseradish peroxidase conjugate is added to all wells. The conjugate binds to IgG / IgM / IgA antibodies respectively on the solid phase antigen to form a stable sandwich. After a 30 minutes incubation at room temperature the excess conjugate is removed by washing all wells with washing buffer.

#### 3.3 Substrate reaction and stopping

The TMB substrate is dispensed into each well and the peroxidase enzyme/substrate reaction forms a stable blue chromogen. The reaction and subsequently the colour development is stopped after 20 minutes incubation at room temperature by adding 0.5 N H<sub>2</sub>SO<sub>4</sub> to the wells. The change in pH also causes the chromogen to change colour from blue to yellow.

#### 3.4 Reading and interpretation

The intensity of the colour is read in a microtiter plate reader at 450 nm (recommended reference wavelength for bichromatic measurement: 600 - 690 nm). The intensity of the colour (OD) is directly proportional to the concentration of the specific antibody in the patient sample.

### 4. Kit Contents

The kits contains sufficient reagents for 12 x 8 = 96 determinations. The strips and solutions have to be stored at 4 - 8 °C. The expiry date is mentioned on the labels.

- 12 strips            Microtiter strips                            single strips each with 8 break-apart wells coated with antigen of HSV-1/2
- 1 x                    Frame holder
- 4 x 2 mL            Calibrators 1 - 4                                    human serum containing antibodies against HSV-1/2 (concentrations listed below) diluted in PBS and stabilised with 0.01 % methylisothiazolone and 0.01 % bronitrodioxane as preservatives, ready to use.

		IgG	IgM	IgA
Cal. 1 (negative)	Concentration (U/mL)	1	1	1
Cal. 2 (cut-off)		10	10	10
Cal. 3 (weak positive)		50	50	50
Cal. 4 (positive)		200	100	100

- 1 x 60 mL            Serum diluent                                    PBS/BSA buffer solution, contains < 0.1 % sodium azide as preservative, ready to use
- 1 x 12 mL            Enzyme conjugate solution                    HRP-labelled goat anti-human-IgG /-IgM /-IgA, ready to use
- 1 x 12 mL            TMB substrate                                    3,3',5,5' Tetramethylbenzidine, ready to use
- 1 x 12 mL            Stopping solution                                0.5 N sulfuric acid, ready to use
- 1 x 60 mL            Washing buffer                                    PBS/Tween buffer solution 10x concentrated to be diluted 1:10 prior to use; the concentrate should be warmed up to 37 °C for 15 min to avoid any crystals
- 2 x                    Plate sealers                                      to cover microtiter strips during incubation
- 1 x                    Plastic bag                                         re-sealable for dry storage of non-used strips

### 5. Materials Required but not Provided

- 5 µL-, 100 µL- and 500-µL micro- and multichannel pipets
- Microtiter plate reader with a 450 nm filter (reference filter 600 – 690 nm)
- Microtiter Plate Washer (in case of manual washing: wash bottle)
- Reagent tubes for the serum dilution
- Measuring cylinder
- Distilled water or water of higher quality

## 6. Warning and Precautions

- For *in-vitro* diagnostic use only! Do not ingest or swallow! Laboratory safety precautions should be followed. Do not eat, drink or smoke in the laboratory.
- All sera and plasma or buffers based upon have been tested to HBsAg, HIV and HCV respectively with generally accepted methods and were found negative. Nevertheless, precautions like the use of latex gloves have to be taken.
- Serum and reagent spills have to be wiped off with a disinfecting solution (e.g. sodium hypochlorite 5 %) and have to be disposed of properly.
- All reagents have to be brought to room temperature (18 to 24 °C) before performing the test.
- Before pipetting all reagents should be mixed thoroughly by gentle tilting or swinging. Vigorous shaking with formation of foam should be avoided.
- It is important to pipet with constant intervals so that all the wells of the microtiter plate have the same conditions.
- When removing reagents out of the bottles care has to be taken that the stoppers are not contaminated. Further a possible mix-up has to be avoided. The content of the bottles is usually sensitive to oxidation so that they should be opened only for a short time.
- In order to avoid a carry-over or a cross-contamination separate disposable pipet tips have to be used.
- No reagents from different kit lots should be used and they should not be mixed with one another.
- All reagents have to be used within shelf life.
- In accordance with a Good Laboratory Practice (GLP) or following ISO 9001 all laboratory devices employed should be regularly checked regarding the accuracy and precision. This refers e.g. to microliter pipets and washing or reading (ELISA Reader) instrumentation.
- The contact of certain reagents especially the stopping solution and the substrate with skin, eye and mucosa has to be avoided because possible irritations and acid burns could arise and there exists a danger of intoxication.

## 7. Storage and Stability

Store all reagents at 4 - 8 °C.

The expiry date of each reagent is printed on the individual labels. Do not use any reagents after the expiry date has been exceeded.

The diluted washing buffer is stable for up to 4 weeks when stored at 4 - 8 °C.

## 8. Specimen Collection and Handling

Both serum or plasma (EDTA, heparin) can be used for the determination. Serum is separated from the blood which is aseptically drawn by venipuncture after clotting and centrifugation. The serum or plasma samples can be stored at 4 - 8 °C for up to 48 hours. They should be kept at -20°C for a longer storage. The samples should not be frozen and thawed repeatedly. Lipemic, hemolytic or bacterially contaminated samples can cause false positive or false negative results.

**Patient sera** must be prediluted **1:101** in **serum diluent** (e.g. 5 µL serum + 500 µL serum diluent) prior to testing.

Samples containing concentrations higher than the highest calibrator have to be diluted further with serum diluent.

In case of interference with rheumatic factors, serum preabsorption with Rf absorbant (MASTSORB Order Code: 651003) is recommended. Do **not** absorb the calibrators.

## 9. Assay Procedure

### 9.1. Preparation of Reagents

**Allow all kit components and specimens to reach room temperature (RT, 18 – 24 °C) prior to use and mix well.**

**Washing buffer:** Dissolve any crystals which may be in the bottle by warming to 37 °C and then mix well.

Dilute the concentrated **washing buffer 1:10 with distilled water** (e.g. 60 mL buffer concentrate + 540 mL distilled water). Mix thoroughly.

- Strictly follow the instructions for reliable test performance. Any changes or modifications are within the responsibility of the user.
- All reagents and samples must be brought to room temperature before use, but should not be left at this temperature for longer than necessary.
- A standard curve should be established with each assay.
- Put the unused microtiter strips back in the plastic bag and store them dry at 4 - 8 °C.

### 9.2. Assay Steps

Prepare a sufficient amount of microtiter wells for calibrators, controls and samples.

**Note: Other incubation conditions might be possible. In case of modifications of the recommended test procedure (e.g. incubation temperature 37 °C instead of RT) the user has to validate assay performance.**

1. Pipette **100 µL** each of the **diluted (1:101) samples** and the **ready to use calibrators** into the appropriate wells.
2. Cover plate with the enclosed plate sealing foil and **incubate** at room temperature for **60 minutes**.
3. Discard the contents of the microwells and wash **3 times** with **300 µL** of **diluted washing buffer**. Afterwards remove residues of the washing solution by gentle tapping of the microtiter plate on a paper towel.
4. Pipette **100 µL** of **enzyme conjugate** solution into each well.
5. Cover plate with plate sealing foil and **incubate** for **30 minutes** at room temperature.
6. Discard the contents of the microwells and wash **3 times** with **300 µL** of **diluted washing buffer**. Afterwards remove residues of the washing solution by gentle tapping of the microtiter plate on a paper towel.
7. Dispense **100 µL** of **TMB substrate** into each well.

8. Cover plate with the plate sealing foil and **incubate** for **20 minutes** in the dark (e.g. drawer) at room temperature.
9. Add **100 µL** of **stopping solution** to each well.
10. After thorough mixing and wiping the bottom of the plate, **read the optical density at 450 nm** and calculate the results. Blank against air. A bichromatic measurement using a reference wavelength of 600 - 690 nm is recommended.

**The developed colour is stable for at least 60 minutes. Read optical densities during this time.**

## 10. Results and Interpretation

### Example

	OD 450 nm	corrected OD	Mean OD Value
Blank	0.020		
Calibrator 1 (negative)	0.150 / 0.144	0.130 / 0.124	0.127
Calibrator 2 (cut-off)	0.530 / 0.570	0.510 / 0.550	0.530
Calibrator 3 (weak positive)	1.100 / 1.132	1.080 / 1.112	1.096
Calibrator 4 (positive)	1.500 / 1.590	1.480 / 1.570	1.525

The table above should be considered as an example which was achieved under arbitrary temperature and environmental conditions. These data do NOT describe **reference values** which have to be found in other laboratories in the same way!

### 10.1. Qualitative Calculation

The calculated OD values for patient sera as mentioned above are compared with the value for the cut-off calibrator. If the value of the sample is higher, then it should be read as positive.

A value below the cut-off calibrator should be read as negative. It seems reasonable to define a range of  $\pm 20\%$  around the value of the cut-off as a grey zone. It is recommended to repeat results laying within the grey zone using the same serum or a new sample of the same patient, taken after 2 - 4 weeks. Both samples should be measured in parallel in the same run.

The positive calibrator must show an absorption value at least double the value received by the cut-off calibrator.

### 10.2 Quantitative Calculation

The ready to use calibrators of the HSV-1/2 antibody ELISAs are defined and values expressed are in arbitrary units (U/mL). This gives access to an exact and reproducible quantification and in consequence patient antibody titer monitoring is possible. Concentration values for calibrators are printed on the labels of the vials.

A standard curve is plotted by entering the mean absorbance value of the calibrators on the Y-axis and the corresponding concentration on the X-axis using graph paper. The concentration of the patient samples can then be read directly from the graph.

The calculation of the result can be performed using a computer and a suitable software program.

## 11. Assay Performance

The assay characteristics of the MASTAZYME HSV-1/2 IgG / IgM / IgA ELISAs have been established and assessed according to the European IVD directive. Detailed validation data can be provided on special request.

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Laboratoriums-Präparate GmbH
**MAST GROUP Ltd.**