MASTAFLUOR™ Borrelia

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Introduction

Borrelia burgdorferi belongs to the family of spirochetes of which three types have been identified to be human pathogenic: *Borrelia burgdorferi sensu stricto*, *B. afzelii* and *B. garinii*. The endemic areas of *Borrelia* are in Central and Eastern Europe as well as Russia, China and Japan. The illness is transferred via tick bites, in Europe mainly by *Ixodes ricinus*. In the endemic zones, like Southern Germany and Austria, up to 50 % of the ticks are infected.

In the clinical course after an erythema migrans, e.g. with neuroborreliosis which appears at the first stage, also chronic arthritis, encephalitis, meningitis, myositis and hepatitis are observed. Treatment is done via different antibiotics, e.g. doxycyclin, amoxicillin, cefuroxim and penicillin G. A specific immunization is possible with immunoprophylaxis either by a recombinant OspA or by a recombinant polyvalent OspC vaccine.

The laboratory diagnosis is performed by the detection of antibodies in blood and cerebrospinal fluid. Methods employed are: ELISA, immunofluorescence, hemagglutination or Western blot. It could be shown that between the various test methods there appear significant differences in the interpretation so that the most reliable method seems to be the follow-up of the titre development. Western blot serves as a confirmatory test because electrophoretically separated single antigens can be evaluated in their reaction with specific serum antibodies.

If the test shows only a positive IgG result (IgM negative) the immunity of the patient can be concluded. In case of a clearly positive clinical finding, confirmation of the titre increase should be performed with the same assay. There is a strong suspicion of an active disease with a simultaneous positive IgM.

Intended Use

MASTAFLUOR™ *Borrelia* provide an indirect immunofluorescence test system for the detection of specific antibodies to *Borrelia burgdorferi*.

Test Principle

IFA slides are coated with purified antigen. Diluted patient serum or plasma is incubated on the wells of the slide. In the presence of specific antibodies a stable antigen-antibody complex is formed. Unspecific or unbound antibodies are removed in a washing step. Specific immune complexes are then detected by either a FITC-conjugated anti-human IgG or IgM antibody. Again unspecific or non-bound antibodies are removed in a washing step. Results could be read visually using an immunofluorescence microscope.

Kit Contents

| Content | 50 Tests | Specifications |
|------------------|--------------|---|
| Slides | 10 x 5 Wells | Coated with purified <i>Borrelia burgdorferi</i> |
| Positive control | 0.5 mL | IgG or IgM positive serum, ready to use, contains < 0.1 % NaN ₃ |
| Negative control | 0.5 mL | IgG and IgM negative serum, ready to use, contains < 0.1 % NaN ₃ |
| FITC-conjugate | 2 mL | Anti-human IgG (γ-chain) FITC-conjugate or anti-human IgM (μ-chain) FITC-conjugate; ready to use, contains < 0.1 % NaN ₃ |
| Serum diluent | 1 x 10 mL | Serum diluent, 5 x conc., contains < 0.1 % NaN ₃ |
| Mounting medium | 3 mL | Ready to use, contains < 0.1 % NaN ₃ |
| PBS | 2 Sachets | 1 sachet to be diluted in 1 L of distilled water, pH 7.2 ± 0.2 |

Materials Required but not Provided

1. Sterile test tubes
2. Micropipettes and tips
3. Staining dish or Coplin jar
4. Moist chamber
5. Volumetric flask for PBS
6. Distilled water or water of higher quality
7. Forceps
8. Transmitted or epi-fluorescent microscope with a 490 nm excitation filter and a 510 nm barrier filter
9. Cover slips
10. Wash bottle

Stability and Storage

MASTAFLUOR™ *Borrelia* is stable until end of shelf life as indicated on kit label. The kit should be stored at 2–8 °C.

Reconstituted PBS powder should be stored at 2–8 °C for up to 30 days.

Serum and plasma specimens may be stored at 2–8 °C for up to 3 days prior to use or for longer storage at -20 °C. Repeated freezing and thawing of specimens should be avoided.

Warnings and Precautions

1. The reagents supplied in this kit are for *in vitro* diagnostic use only.
2. Read instructions carefully before conducting the assay. Do not modify the procedure.
3. Do not use kit beyond the expiry date.
4. Wear appropriate protective clothing and use appropriate facilities.
5. Do not mouth pipette.
6. Use disposable plasticware where ever possible. Re-usable glassware should be washed thoroughly and rinsed free of detergents before use.
7. Sera provided have been tested for antibodies to HIV and for HBsAg and found to be negative. However, they should be treated as potentially infectious materials and capable of transmitting diseases. No guarantee is given that the sera are free of infection or microbial contamination.
8. Use distilled water or water of higher quality.
9. Do not mix reagents between different lots, as reagents have been calibrated for each batch.
10. Do not cross-contaminate reagents or interchange caps on bottles. Use a separate pipette or pipette tips for each sample and reagent.
11. Do not allow wells to dry out during the assay procedure.
12. Do not expose slides to intense sunlight or similar adverse conditions while incubating.
13. Contaminated plastic ware should be disposed and incinerated. Contaminated glassware should either be autoclaved at 121 °C or decontaminated with a solution of 2.5 % sodium hypochlorite. All liquid waste should be decontaminated appropriately e.g. using sodium hypochlorite.
14. Do not expose the FITC-conjugate at any stage to strong sunlight, UV or fluorescent light. Keep in a dark place whenever possible.
15. Microbial contaminated serum samples should not be used.
16. In the event that haemolysed or lipemic serum must be used, heat inactive at 56 °C before use.
17. Sodium azide is used as a preservative as marked. It may be toxic if ingested. Sodium azide may react with lead and copper plumbing to form highly explosive salts. Always dispose by flushing to drain with plenty of water.

Test Procedure

1. Allow all materials to reach room temperature (at least 20 °C) prior to use.
2. Reconstitute one sachet of PBS powder with 1 litre of distilled water or water of higher quality.
3. Dilute concentrated serum diluent 1:5 with distilled water e.g. 20 mL diluent + 80 mL distilled water.
Dilution of patient samples:
For IgG antibodies:
1:160 with serum diluent
e.g. 5 µL serum + 795 µL serum diluent
For IgM antibodies:
1:40 with serum diluent
e.g. 5 µL serum + 195 µL serum diluent
For IgM determination a preabsorption with Rf absorbent (MASTSORB™ Order Code: 651003) is recommended to reduce interferences e.g. with rheumatic factors. **Do not absorb the controls.**
4. Controls are supplied ready to use and must not be diluted.
5. Carefully remove the required number of slides from their sachets and mark accordingly. Hold the slide at the edge and do not touch the wells. Place the slides in a moist chamber.
6. Apply 20–25 µL of diluted specimens and ready to use controls to respective wells on the slides according to a prepared work format. Ensure that all the wells are covered and that serum does not escape from the wells.
Direct contact of a pipette tip with the slide surface may result in damage to the antigen substrate and should be avoided.
7. Cover the moist chamber and incubate as follows:
IgG assay: room temperature for 30 min.
IgM assay: room temperature for 45 min.
8. After incubation wash slides carefully with PBS from a wash bottle, taking care not to direct the jet onto the test wells. This may be done by directing the jet of PBS along the centre of the slide, tilting the slide first towards wells 1–5 and then towards 6–10.
9. Immerse the slides in a staining dish or Coplin jar containing PBS and wash for 15 min with a change of PBS after 5 min.
10. Remove slides one at a time from the staining jar and drain off any excess PBS. Using a blotter, dry the area outside wells. Wipe off any excess PBS from the back of the slide. Do not touch the well surface.
11. Immediately transfer the slide to a moist chamber and add 20–25 µL of FITC-conjugate (either IgG or IgM) to each well. Ensure that all the wells are covered by the conjugate.
12. Incubate slides at room temperature in a covered moist chamber in the dark for 30 minutes (same incubation for IgG and IgM).
13. Switch on microscope 15 minutes prior to use.
14. After incubation wash the slides with PBS as in steps 9 and 10.

15. Apply one small drop of Mounting Medium to each well on a slide. Place a cover slip on the slide and read the results immediately. Take care not to trap air pockets under the cover slip. Excess Mounting Medium on a slide may result in high background fluorescence due to light scattering. Remove excess Mounting Medium with a paper towel avoiding any direct movement of the cover slip.
16. Examine under a fluorescence microscope at a total magnification of 400x to 800x.

Interpretation of Results

Test Validation

Examine the control samples before reading patient samples. The controls must give the following patterns to be valid.

Specimen sera are assessed by way of comparison with the controls included in the test.

| Sample | Fluorescence Intensity |
|----------------------|------------------------|
| Positive control IgG | 3 + to 4 + |
| Positive control IgM | 2 + to 3 + |
| Negative control | negative |

Interpretation of Specimen Results

IgG positive specimen should show a homogenous green fluorescence of all active Borrelia. Less active Borrelia may also be present with a faint green colour and these should not be taken into consideration.

IgG negative specimen should show no fluorescence; a green staining without fluorescence or weak fluorescence (mainly at the edge of the test-field) may be visible.

IgM positive specimen should show a homogenous green fluorescence of all active Borrelia. Less active Borrelia may also be present with a faint green colour and these should not be taken into consideration.

IgM negative specimen should show no fluorescence; a green staining without fluorescence or weak fluorescence (mainly at the edge of the test-field) may be visible.

IgG-Ab \geq 1:160 and IgM-Ab \geq 1:40

Previous infection probable, possible reinfection or beginning infection.

IgG-Ab \geq 1:320 and IgM-Ab \geq 1:80

Acute or recent reinfection probable, reactivation possible.

Limitations / Cross-reactions

1. Cross-reactivity with spirochetes, particularly *T. pallidum* can result in false positives. In case of positive and borderline results it is recommended to exclude infections with cross reacting organism.
2. Possible cross-reactivity may be avoided by using spirochete absorbent.
3. In pre-absorbed samples detection of low Borrelia antibody titres may be not detected.
4. MASTAFLUOR™ Borrelia has been demonstrated to be highly sensitive and specific, however, results should be considered along with all serological tests, clinical history and other aspects of patient management to be considered diagnostically significant.

Performance

Data in the table were obtained by comparing MASTAFLUOR™ Borrelia results with data of other assays based on a different method (e.g. ELISA):

| | IgG | IgM |
|-------------|---------|--------|
| Sensitivity | 95.2 % | 57.1 % |
| Specificity | 98.75 % | 97.8 % |

Inter-assay and intra-assay variances are not detectable in respect to control fluorescence intensity.

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